

## Highlights

3 recording modes: auto-trigger (detection of tachycardia, bradycardia, atrial fibrillation and pauses), time based and manually triggered

Selectable pre- and post-storage durations

Operation up to 14 days with one battery

Snap-on standard gel electrodes – no need for ECG cables

Mobile ECG-transmission via PhysioGate® App

Ingress protection against dust and splashing water (IP 64)

# CardioMem® CM 100 XT

## External ECG-Loop Recorder incl. PhysioGate® App

The CardioMem® CM 100 XT is a mobile ECG loop recorder designed to acquire and record ECG stripes in three ways: upon manual activation by the patient, upon automatic detection of arrhythmic events (tachycardia, bradycardia, atrial fibrillation and pauses), or at predefined intervals via an integrated timer. The physician can configure the detection parameters for each patient individually. This versatility enables the physician to diagnose the patient effectively.


### Easy to use – versatile

Standard adhesive electrodes ensure maximum comfort for the patient without the burden of cumbersome lead wires. Application of the device can be easily performed directly by the patient by simply removing the electrode or patch protection foil and pressing the device onto the chest in close proximity to the heart. The auto-trigger algorithm for detecting tachycardia, bradycardia, atrial fibrillation and pauses as well as the interval timer for automatic signal acquisition at predetermined intervals make the CardioMem® CM 100 XT the optimal device for detecting and recording rare and asymptomatic cardiac arrhythmias. The ability to specifically select the types of events that should be detected and the appropriate pre- and post-times for these events further enhance the device's versatility.

### ECG quickly available via PhysioGate® App

The data are automatically transmitted to the physician's office or a clinic via the PhysioGate® App, thus ensuring that the ECG events are readily available while minimizing the length of time the patient needs to wear the device.

The CardioMem® CM 100 XT in connection with the PhysioGate® App aids in performing a fast diagnosis and early initiation of treatment, thereby contributing to increase the patient's quality of life. The robust design ensures that it can be used over and over again by a large number of patients.




# CardioMem® CM 100 XT

## Technical Specifications

<b>General</b>	
Dimensions	76 mm x 89 mm x 14 mm
Weight	approx. 39 g (including battery)
Battery type	3 V Lithium, Renata CR2477N
Operation time	7 ... 14 days
Operation	Event button
Connections	Snap receptacles for standard adhesive electrodes
Others	Port for custom USB cable
<b>Parameters</b>	
ECG leads	2 channels, 3 electrodes
Frequency range	0,05 ... 70 Hz
Open-Lead recognition	Yes
Time-triggered recording	1–24 h (interval 1 h)
Pauses	2–4 s (interval 0,5 s)
Atrial fibrillation	On / Off
Bradycardia	30–60 1/min
Tachycardia	100–240 1/min
Pre-time (selectable)	25 / 30 / 45 / 60 s
Post-time (selectable)	15 / 30 / 60 s
<b>Storage</b>	
Digitale signal processing	512 Hz / 24 Bit
<b>Data transmission</b>	
Transmission method	PhysioGate® App
<b>Classification</b>	
Product classification	Ila according to (EU) 2017/745 (MDR)
Classification of the applied part	Type BF (Body Floating)
Ingress protection	IP 64
UMDNS-Code	12-601

<b>Operation conditions</b>	
Temperature	5 ... 45 °C
Relative humidity	0 ... 93 %, not condensing
<b>Storage and transport conditions</b>	
Temperature	-20 ... 60 °C
Relative humidity	0 ... 93 %, not condensing
<b>Radio Power</b>	
Radio technology	Bluetooth Low Energy
Frequency	2 402 MHz ... 2 480 MHz (2,4 GHz ISM Band)
Modulation	Gaussian frequency shift keying (GFSK)
Radiation	0 dBm = 1 mW



Subject to change



Gefördert durch: Bundesministerium für Wirtschaft und Energie aufgrund eines Beschlusses des Deutschen Bundestages

Manufacturer



GETEMED  
 Medizin- und Informationstechnik AG  
 Oderstr. 77 / 14513 Teltow  
 Tel. +49 3328 3942-0  
 Fax +49 3328 3942-99  
 info@getemed.de / www.getemed.de